

# MONARCH™

# 21

## A Q U A N T U M L E A P

MONARCH 21 takes penile plethysmography (PPG) into the 21st century with enhanced flexibility, portability and convenience.



**Flexibility:** MONARCH™ 21 Runs Standardized and Non-Standardized Assessments!

MONARCH™ 21 runs standardized assessments and non-standardized PPG assessments. In addition, MONARCH™ 21 standardized assessments can include “specialty segments” randomly integrated into a standardized stimulus presentation... all defined by the clinician. The non-standardized assessment option enables the clinician to run virtually any stimulus, providing more flexibility for research, assessment and conditioning.

**Portability: Palm-sized DRD!**

The MONARCH™ 21 Data Recording Device (DRD) is small enough to fit in the palm of your hand and it only has three wires! That’s right...only THREE. Additionally, digital stimulus delivery is controlled and delivered in software. If you use compatible video goggles you won’t even need a TV monitor. MONARCH™ 21 is a truly simplified and portable system!!

**Convenience & Versatility: MONARCH™ 21 Software Runs in Windows XP! MONARCH™ 21 Male Assessment Software is extremely intuitive**

and easy to use with a clinician-designed user interface. Real-world PPG administration and application experience has gone into the design of MONARCH™ 21’s intuitive and easy-to-understand software and hardware.

MONARCH™ 21 is the latest in a long line of PPG technology and treatment breakthroughs: the first standardized PPG system, the first research-based and appropriate PPG stimuli (adults, adolescents and mentally-handicapped offenders), effective anti-faking features, and the lowest non-responder rate in the field. MONARCH™ 21 maintains these gains, while adding to them the technological advantages described

above. Contact us to learn how you can put “the System that Works!” to work for you as an indispensable part of the treatment and assessment services you offer.



MONARCH™ 21 palm-sized Data Recording Device (DRD)

# MONARCH™ 21 System includes:

MONARCH™ 21 Data Recording Device. Three wires collect three simultaneous physiological data streams... Penile, Respiration and GSR... for a more detailed and reliable PPG. No more knobs, dials or buttons. Set-up and calibration are a snap with on-screen software controls.

Easy-to-use MONARCH™ 21 PPG Software. The world's only standardized PPG with MONARCH™ 21 DVD stimulus sets. MONARCH™ 21 is also user-definable and can be used with customized and specialized stimulus for investigation and research. MONARCH™ 21 is configured on a state-of-the-art laptop PC.

MONARCH™ 21 Intercom Link connects technician and client with clear voice and audio communication.

Three sizes of Barlow and Indium Gallium gauges. Includes precision calibration cone set.



Full duplex technician and client head sets with integrated boom mics.

MONARCH™ 21's Respiration, GSR sensors and Event Marker Button enhance client participation while collecting valuable test behavior data.



Optional video goggles maximize clarity of MONARCH™ 21 DVD visual stimulus (call for price and availability).



**Your choice...** of Adult, Adolescent or DDMR non-pornographic DVD stimulus sets. MONARCH™ 21 stimulus sets are the only clinically validated, research-proven stimulus available anywhere. MONARCH™ 21 software randomizes each PPG session using MONARCH™ 21 DVD stimuli. Specialty segments are included and can be easily added to test client-specific issues in assessment or treatment. MONARCH™ 21 stimuli are digitally mastered in DVD format for maximum audio and visual fidelity.

**Four days of training . . .** taught by PPG experts with real-world clinical experience. MONARCH™ 21 purchasers are given two days of technician and two days of clinical training at Behavioral Technology in Salt Lake City. Instruction by clinical experts means you'll be interpreting meaningful PPG assessments sooner.

The bottom line... *Everything you need to do the best PPG.*